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
The screenshot displays the AutoTEM Cryo software interface, which is organized into several panels. On the left is the 'TEMPLATE LIST' panel, showing a search bar and a list of templates: 'Yeast default template Cryo' (highlighted with a red box), 'auto_sleep Cryo', and 'Cell Lamella template Cryo'. The main area is divided into three vertical sections: 'Preparation', 'Milling', and 'Thinning'. The 'Preparation' section includes settings for 'Ion HFW Oversize' (30.0 μm), 'Sample Pretilt' (45.0 $^\circ$), and 'Milling angle' (Not Available). The 'Milling' section includes 'Size' (5.0 μm x 3.0 μm), 'Correction Factor' (0.50), and a list of steps with toggle switches: DELAY, ELECTRON REFERENCE DEFINITION, REFERENCE DEFINITION, STRESS RELIEF CUTS, REFERENCE REDEFINITION 1, ROUGH MILLING, ROUGH MILLING - ELECTRON IMAGE, REFERENCE REDEFINITION 2, MEDIUM MILLING, MEDIUM MILLING - ELECTRON IMAGE, FINE MILLING, FINE MILLING - ELECTRON IMAGE, FINER MILLING, and FINER MILLING - ELECTRON IMAGE. The 'Thinning' section includes 'Final Thickness' (250.0 nm), 'Enable Windows' (disabled), and a list of steps with toggle switches: DELAY, POLISHING 1, POLISHING 1 - ELECTRON IMAGE, POLISHING 2, and POLISHING 2 - ELECTRON IMAGE. A red box highlights the 'EUCENTRIC TILT' toggle switch in the 'Preparation' section, which is currently turned on. A red text overlay in the bottom left corner of the 'Preparation' section reads: '采用倾转样品台寻找共焦点，耗时长、样品台损耗大而且精度不理想'. A yellow box highlights the 'Thinning' section, with a yellow text overlay in the center reading: '非必要'. The bottom status bar shows the version '2.2.0 (core: 10.0.3)'.

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

AutoScript

AutoScript

AutoScript is a software suite that allows to *perform microscopy analysis and related data processing tasks* in an Python 3 language environment.

The suite contains Python 3 distribution with an interpreter and a rich set of widely used packages like `opencv-python`, `numpy` or `scipy`. Most importantly, it also contains a set of special packages related to microscopy analysis, particularly `autoscript_sdb_microscope_client` with **SDB microscope control features** and `autoscript_toolkit` with **domain-specific data processing capabilities** like fine-tuned computer vision algorithms.

This is the **reference manual** for both forementioned *essential packages*. Please explore the following sections to find detailed information:

[SDB microscope client](#)

[Toolkit](#)

A **custom task library** with a minimum number of dependencies.

Numpy & SciPy

- Linear algebra & Scientific computing

OpenCV

- Computer Vision Library for real-time image processing

QT4

- Design of simple dialogs to move away from the command line with minimal effort.

Matplotlib

- Plotting images and charts

Libraries often used but not in custom library.

Pandas

- Database / Data analysis

SciKit Image

- Image processing / Machine learning

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Name	Date modified	Type
7degree	09/02/2023 15:23	File folder
8degree	09/02/2023 15:21	File folder
9degree	09/02/2023 15:19	File folder
10degree	03/02/2023 10:32	File folder
11degree	09/02/2023 15:22	File folder
12degree	09/02/2023 15:24	File folder
13degree	09/02/2023 15:25	File folder
14degree	09/02/2023 15:26	File folder
15degree	09/02/2023 15:27	File folder
16degree	09/02/2023 15:28	File folder
17degree	09/02/2023 15:29	File folder
18degree	09/02/2023 15:29	File folder
19degree	09/02/2023 15:30	File folder
20degree	09/02/2023 15:30	File folder

Name	Date modified	Type	Size
default_json	09/02/2023 15:23	File	1 KB

user_folder

D:/SharedData/XiaoJunHuang/auto_script/auto_slice_example

set reference

take map image

load map

name, stage_x, stage_y, stage_z, milling_angle, thinned

point_001 3.49 2.89 6.92 10.0 N

add point

update point

E correct

delete point

goto point

update pattern

settings for taking image

electron

Magnification

350x

high tension

3.0 kV

beam current

0.10 nA

dwelt time

1 us

resolution

1536x1024

ion

2000x

30.0 kV

0.10 nA

100 ns

1536x1024

auto milling current position

auto milling all points

stop

settings for milling with ion beam

name, position, type, high tension, beam current, milling angle, x, y, z

add pattern

change pattern

delete pattern

include side

show pattern

load pattern

manual thin toolset

tilt to milling angle 1

15.3

show select pattern

ccs side pattern

tilt to milling angle 2

14.7

clear pattern

+ pattern

take E image

save E image


save I image

take I image

pause

stop

start milling

user_folder D:/SharedData/XiaoJunHuang/auto_script/auto_slice_example

set reference take map image load map

name, stage_x, stage_y, stage_z, milling_angle, thinned

point_001 3.49 2.89 6.92 10.0 N

add point update point E correct delete point goto point update pattern

settings for taking image

	Magnification	high tension	beam current	dwel time	resolution
electron	350x	3.0 kV	0.10 nA	1 us	1536x1024
ion	2000x	30.0 kV	0.10 nA	100 ns	1536x1024

=====

auto milling current position auto milling all points stop

=====settings for milling with ion beam=====

name, position, type, high tension, beam current, milling angle, x, y, z

add pattern change pattern delete pattern ☒ include side show pattern load pattern



=====manual thin toolset=====

tilt to milling angle 1	15.3	show select pattern	ccs side pattern
tilt to milling angle 2	14.7	clear pattern	+ pattern
take E image	save E image	save I image	take I image
	pause	stop	
start milling			

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






user_folder
D:/SharedData/XiaoJunHuang/auto_script/auto_slice_example

set reference
take map image
load map

name, stage_x, stage_y, stage_z, milling_angle, thinned

point_001	3.49	2.94	7.00	10.0	N
point_002	3.61	2.92	6.99	10.0	Y
point_003	3.49	2.94	7.00	10.0	N

add point
update point
E correct
delete point
goto point
update pattern

settings for taking image

	Magnification	high tension	beam current	dwel time	resolution
electron	350x	3.0 kV	0.10 nA	1 us	1536x1024
ion	2000x	30.0 kV	0.10 nA	100 ns	1536x1024

=====

auto milling current position
auto milling all points
stop

=====settings for milling with ion beam=====

name, position, type, high tension, beam current, milling angle, x, y, z

pattern_001	bottom	rectangle	30.0 kV	0.3 nA	10.0degree	13um	5.3um	0.35um
pattern_002	Top	rectangle	30.0 kV	0.3 nA	10.0degree	13um	5.0um	0.35um
pattern_003	bottom	rectangle	30.0 kV	0.10 nA	9.5degree	11um	1.5um	0.35um
pattern_004	Top	rectangle	30.0 kV	0.10 nA	10.5degree	11um	1.5um	0.35um

add pattern
change pattern
delete pattern
☒ include side
show pattern
load pattern

=====manual thin toolset=====

tilt to milling angle 1	15.3	show select pattern	ccs side pattern
tilt to milling angle 2	14.7	clear pattern	+ pattern
take E image	save E image	save I image	take I image
pause	stop		
start milling			

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